

# **JEDNORÁZOVÝ RYCHLOTEST** **na detekci COVID-19**

Souprava k detekci hrotového glykoproteinu  
nového koronaviru (ligand-receptorová  
kompetitivní chromatografie)



Obsah testovací sady Newgene

## VÝHODY ANTIGENNÍCH TESTŮ

Použití vzorku sputa (hlen)

Rychlá detekce - výsledek do 15-ti minut

Vysoká přesnost - klinická účinnost 95.8%

Snádné a nebolestivé použití

## ZOBRAZENÍ VÝSLEDKŮ

Výsledek se zobrazí do 15-ti minut



Pozitivní




Negativní



Neplatný

## POSTUP PŘI TESTOVÁNÍ

Vzorek ze sputa



Pomocí vatové tyčinky naberte 10 až 50 mg vzorků sputa (odpovídá velikosti zápalkové hlavičky).

Otevřete krytku extrakční zkumavky na vzorek a ulomte špičku stěrové tyčinky do zkumavky. Zavřete jednorázovou extrakční zkumavku na vzorek a protřepáváním nechte zcela promíchat.

Vytáhněte z balicího sáčku testovací kazetu, uložte ji na stůl, odlomte vyčnívající část odběrové zkumavky a vertikálně přidejte 3 kapky vzorku do otvoru na vzorek.



EAN	8456713601546
Počet kusů v balení	1x testovací sada
Doba dodání	Do 14-ti dní



Medicines & Healthcare products  
Regulatory Agency



**Our Ref: IVD001178**

Dr Edward Wang  
Wellkang Ltd  
16 Castle Street  
Dover  
Kent  
CT16 1PW  
United Kingdom

**MHRA**

10 South Colonnade  
Canary Wharf  
London  
E14 4PU  
United Kingdom

**18 May 2020**

Dear Dr Wang

**IN VITRO DIAGNOSTIC MEDICAL DEVICES REGULATIONS 2002: REGULATION 44**  
**Registration of manufacturers of *In-Vitro Diagnostic* Medical Devices**  
**and devices for Performance Evaluation**

Thank you for informing the Competent Authority of the details of **Manufacturers Name:- New Gene (Hangzhou) Bioengineering Co., Ltd.** located at **Manufacturers Address:- Room 1606, Floor 16, Building 5, 688 Bin'an Road, Changhe Street, Binjiang District, Hangzhou City, Zhejiang, China 310000** for whom you are acting as the authorised representative and for supplying the medical device information.

**Your registration has been recorded based on your declaration that you have determined that the device(s) fall within the definition of "in vitro diagnostic medical device", and that you have classified it/them correctly considering the intended purpose(s) and mode(s) of action. In accepting your registration, I should make clear that the Competent Authority does not examine each individual notification and therefore cannot and does not necessarily endorse these determinations.**

**Please note this letter does not represent any form of accreditation, certification or approval by the UK Competent Authority.**

**If you stop placing devices on the market or if you are not complying with the Regulations you should inform us as required by the Regulations. You should be aware that it is an offence to place on the market CE marked devices that do not comply with the Regulations.**

The information you provided has been recorded against the reference number shown at the top of this letter, which we ask you to quote in all future correspondence and communications.

**Please inform us of any of the following changes;**

- the company information
- additional generic groups of devices or, for Annex II or Self-Test devices, additional devices
- discontinuation of a generic group of devices or, for Annex II or Self-Test devices, discontinuation of devices

You should submit your change of registration via DORS with the required statutory fee, which should be accompanied with the information when it is supplied, (the fee is payable for each record notified, and you may place multiple changes on one record).

Thank you for registering the following generic groups of devices

- 1. Part 5: IVDs which are not Annex II and not self-test devices**
- 2.**



Medicines & Healthcare products  
Regulatory Agency



3. *For reagents, reagent products, calibration and control materials:*
4. *group by common technological characteristics and/or analytes*
- 5.
6. *New products:*
7. *None*
- 8.
9. *For performance evaluation:*
10. *None*
- 11.
12. *Neither:*
13. *Coronavirus*
14. *Multiple Drugs of Abuse/Toxicology Rapid Tests*
- 15.
- 16.
17. *For other IVDs, group by appropriate indications*
- 18.
19. *New products:*
20. *None*
- 21.
22. *For performance evaluation:*
23. *None*
- 24.
25. *Neither:*
26. *None*
- 27.
- 28.
29. *Part 6: IVDs which are Annex II or self-test devices*
- 30.
31. *For reagents, reagent products, calibration and control materials:*
32. *group by common technological characteristics and/or analytes*
- 33.
34. *New products:*
35. *None*
- 36.
37. *For performance evaluation:*
38. *None*
- 39.
40. *Neither:*
41. *None*
- 42.
- 43.
44. *For other IVDs, group by appropriate indications*
- 45.
46. *New products:*
47. *None*
- 48.
49. *For performance evaluation:*
50. *None*
- 51.
52. *Neither:*
53. *None*
- 54.



Medicines & Healthcare products  
Regulatory Agency

If you have any queries regarding your registration, please do not hesitate to contact us.



Yours sincerely



[Malcolm Ridgway](#)

Data Integrity Support Officer





**EC Declaration of Conformity**



according to the Directive 98/79/EC  
(applicable to IVD Devices of NOT Annex II and NOT self-test)

**Manufacturer:** **New Gene (Hangzhou) Bioengineering Co., Ltd.**

**Address:** Room 1606, Floor 16, Building 5, 688 Bin'an Road, Changhe Street, Binjiang District, Hangzhou City, Zhejiang Province, P. R. China

**EC Representative:** Wellkang Ltd  
16 Castle Street, Dover, CT16 1PW, England, UK

**We, the manufacturer, declare under our sole responsibility that**

<b>the medical device(s)</b>	<b>Product Name</b>	Novel Coronavirus Spike Glycoprotein Detection Kit (Ligand-Receptor Competitive Chromatography)
	Type/model, identification of product allowing traceability (Where applicable)	COVID-19-NG04
<b>of Category</b>	: <b>Common/Others IVD</b> (Devices of NOT <b>Annex II</b> and <b>NOT self-test</b> )	

**is/are in conformity with the relevant provisions and requirements of Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices.**

Applied harmonised standards, national standards or other normative documents	EN 23640:2015	EN 13640:2002
	EN 980:2016	EN 13641:2002
	EN ISO 14971:2019	EN ISO 18113-1 2011
	EN 13612:2002	EN ISO 18113-4 2011

Conformity assessment procedure	<b>Module A (EC Declaration of Conformity) (Annex III, except point 6)</b>
Notified Body (name & number)	<b>NOT applicable</b>

Signed on: 7 May 2020. Place: Hangzhou City, Zhejiang Province, P. R. China

Signature (on behalf of the manufacturer)

*Mingfu Li*

Name of authorized signatory: Mingfu Li  
Position held in the company: General Manager  
Company Seal/Stamp:



By Royal Charter

# Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that: **New Gene (Hangzhou) Bioengineering Co., Ltd.**  
Room 1606, 16th Floor, No.5 Building  
688 Bin'an Road  
Binjiang District  
Hangzhou  
Zhejiang  
310052  
China

诺迦（杭州）生物工程有限公司  
中国  
浙江省  
杭州市  
滨江区  
长河街道滨安路688号  
5幢16层1606室  
邮编：310052

Holds Certificate No: **MD 729179**

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

Design and Development, Manufacture and Distribution of In-vitro Diagnostic Rapid Test Kit of Drug Abuse, Manufacture and Distribution of In-vitro Diagnostic Rapid Test Kit of Infectious Diseases.

药物滥用体外诊断快速检测试剂盒的设计，开发，制造和销售，传染病体外诊断快速检测试剂盒的制造和销售。



For and on behalf of BSI:

**Gary E Slack, Senior Vice President - Medical Devices**

Original Registration Date: 2020-07-27

Latest Revision Date: 2020-07-27

Effective Date: 2020-07-27

Expiry Date: 2023-07-26

Page: 1 of 1



...making excellence a habit.™